

SNX5 antibody [N2C2], Internal

Cat. No. GTX118542

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human SNX5. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX118542 WB Image

Sample (30 ug of whole cell lysate)


A: A431

B: HeLa

C: HepG2

10% SDS PAGE

GTX118542 diluted at 1:3000




GTX118542 ICC/IF Image

SNX5 antibody [N2C2], Internal detects SNX5 protein at cytoplasm by immunofluorescent analysis.

Sample: A431 cells were fixed in -20°C 100% MeOH for 5 min.

Green: SNX5 protein stained by SNX5 antibody [N2C2], Internal (GTX118542) diluted at 1:500.

Blue: Hoechst 33343 staining.




GTX118542 IHC-P Image

SNX5 antibody [N2C2], Internal detects SNX5 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

SNX5 stained by SNX5 antibody [N2C2], Internal (GTX118542) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX118542 IHC-P Image

SNX5 antibody [N2C2], Internal detects SNX5 protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat tissues.

SNX5 stained by SNX5 antibody [N2C2], Internal (GTX118542) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).