

## Interferon gamma antibody [F1]

Cat. No. GTX11855

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	IP, Immunoassay, Purification
<b>Reactivity</b>	Mouse

Package  
50 µg

## Applications

## Application Note

For IA: Use at an assay dependent dilution. For IP: Use at an assay dependent dilution. For P: Use at an assay dependent dilution. Is unsuitable for Neutralising. Optimal dilutions/concentrations should be determined by the end user.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.1% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant full length protein (Mouse).
<b>Purification</b>	Protein G affinity purified From polyclonal serum
<b>Conjugation</b>	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).