

EFHA1 antibody

Cat. No. GTX118555

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000


Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.55 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human EFHA1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX118555 IHC-P Image

EFHA1 antibody detects EFHA1 protein at cytoplasm in rat colon by immunohistochemical analysis.

Sample: Paraffin-embedded rat colon.

EFHA1 antibody (GTX118555) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX118555 IHC-P Image

EFHA1 antibody detects EFHA1 protein at cytoplasm in mouse kidney by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

EFHA1 antibody (GTX118555) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX118555 IHC-P Image

EFHA1 antibody detects EFHA1 protein at cytoplasm in mouse colon by immunohistochemical analysis.

Sample: Paraffin-embedded mouse colon.

EFHA1 antibody (GTX118555) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX118555 IHC-P Image

EFHA1 antibody detects EFHA1 protein at cytosol on human hepatoma by immunohistochemical analysis.

Sample: Paraffin-embedded hepatoma tissue.

EFHA1 antibody (GTX118555) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX118555 IHC-P Image

EFHA1 antibody detects EFHA1 protein at cytoplasm in rat stomach by immunohistochemical analysis.

Sample: Paraffin-embedded rat stomach.

EFHA1 antibody (GTX118555) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).