

ADI1 antibody

Cat. No. GTX118578

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human ADI1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx118578 ICC/IF Image


Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using ADI1(GTx118578) antibody at 1:200 dilution.


GTx118578 WB Image

Sample (30 ug of whole cell lysate)
A: HepG2
12% SDS PAGE
GTx118578 diluted at 1:2000


GTx118578 IHC-P Image

Immunohistochemical analysis of paraffin-embedded CL1-5 xenograft, using ADI1(GTx118578) antibody at 1:500 dilution.
Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx118578 WB Image

ADI1 antibody detects ADI1 protein by Western blot analysis.
A. 50 µg mouse lung lysate/extract
12 % SDS-PAGE
ADI1 antibody (GTx118578) dilution: 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).