

WBP11 antibody

Cat. No. GTX118654

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
IHC-P	1:100-1:1000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human WBP11. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX118654 IP Image

WBP11 antibody immunoprecipitates WBP11 protein in IP experiments.

IP samples: Jurkat whole cell extract


A. Control with 4 µg of preimmune Rabbit IgG

B. Immunoprecipitation of WBP11 protein by 4 µg WBP11 antibody (GTX118654)


7.5 % SDS-PAGE

The immunoprecipitated WBP11 protein was detected by WBP11 antibody (GTX118654) diluted at 1:1000.

[EasyBlot anti-rabbit IgG (GTX221666-01) was used as a secondary reagent]


GTX118654 WB Image

Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with WBP11 antibody (GTX118654) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX118654 IHC-P Image

WBP11 antibody detects WBP11 protein at nucleus in rat brain by immunohistochemical analysis.

Sample: Paraffin-embedded rat brain.

WBP11 antibody (GTX118654) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX118654 IHC-P Image

WBP11 antibody detects WBP11 protein at nucleus in rat adrenal by immunohistochemical analysis.

Sample: Paraffin-embedded rat adrenal.

WBP11 antibody (GTX118654) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).