

## C3 / C3a antibody [2898]

**Cat. No. GTX11874**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	WB, Immunoassay
<b>Reactivity</b>	Human

**Package**

50 µg

## Applications

**Application Note**

For IA: Use at an assay dependent dilution. For WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

<b>Calculated MW</b>	187 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	The antibody recognizes intact C3, C3a as well as C3a-desArg.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.1% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Unfortunately, this information is considered to be commercially sensitive
<b>Purification</b>	Protein G affinity purified From polyclonal serum
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).