

## GTF3C4 antibody

Cat. No. GTX118874

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human             |

References ( 1 )

Package

100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:3000         |
| ICC/IF             | 1:100-1:1000         |
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.

Calculated MW 92 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 1% BSA, 20% Glycerol                                                                                                                                                                                                  |
| Preservative  | 0.01% Thimerosal                                                                                                                                                                                                           |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Recombinant protein encompassing a sequence within the center region of human GTF3C4. The exact sequence is proprietary.                                                                                                   |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx118874 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using GTF3C4(GTx118874) antibody at 1:500 dilution.



### GTx118874 WB Image

Sample (30 ug of whole cell lysate)  
A: U87-MG  
7.5% SDS PAGE  
GTx118874 diluted at 1:1000



### GTx118874 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human breast cancer, using GTF3C4(GTx118874) antibody at 1:250 dilution.  
Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).