

HRSP12 antibody

Cat. No. GTX119107

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 14 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human HRSP12 Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX119107 ICC/IF Image

HRSP12 antibody detects HRSP12 protein at cytoplasm and nucleus by immunofluorescent analysis.

Sample: HepG2 cells were fixed in ice-cold MeOH for 5 min.

Green: HRSP12 protein stained by HRSP12 antibody (GTX119107) diluted at 1:500.

Blue: Hoechst 33342 staining.



GTX119107 IHC-P Image

HRSP12 antibody detects HRSP12 protein at nucleus on SkHep1 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded SkHep1 xenograft.

HRSP12 antibody (GTX119107) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX119107 WB Image

Sample (30 ug of whole cell lysate)

A: HCT116

15% SDS PAGE

GTX119107 diluted at 1:2000



For full product information, images and publications, please visit our [website](https://www.genetex.com).