

U1SNRNPBP antibody [N2C3]

Cat. No. GTX119546

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 29 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human U1SNRNPBP. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTXT119546 ICC/IF Image

U1SNRNPBP antibody [N2C3] detects U1SNRNPBP protein at nucleus by immunofluorescent analysis.

Sample: U87-MG cells were fixed in 4% paraformaldehyde at RT for 15 min.


Green: U1SNRNPBP protein stained by U1SNRNPBP antibody [N2C3] (GTXT119546) diluted at 1:500.

Blue: Hoechst 33342 staining.


GTXT119546 IHC-P Image

Immunohistochemical analysis of paraffin-embedded Hela xenograft, using U1SNRNPBP(GTXT119546) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTXT119546 WB Image

Sample (50 ug of whole cell lysate)

A: mouse brain

12% SDS PAGE

GTXT119546 diluted at 1:1000


GTXT119546 WB Image

Sample (30 ug of whole cell lysate)

A: U87-MG

12% SDS PAGE

GTXT119546 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).