

## C2orf65 antibody [N2C2], Internal

Cat. No. GTX120030

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB, IHC-P  |
| Reactivity   | Human      |

Package  
100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000-1:10000       |
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.

Calculated MW 59 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | 0.1M Tris, 0.1M Glycine, 20% Glycerol                                                                                                                                                                                      |
| Preservative  | 0.01% Thimerosal                                                                                                                                                                                                           |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Recombinant protein encompassing a sequence within the center region of human C2orf65. The exact sequence is proprietary.                                                                                                  |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX120030 WB Image

Sample (30 ug of whole cell lysate)

A: U87-MG

10% SDS PAGE

GTX120030 diluted at 1:5000



### GTX120030 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human hepatoma, using C2orf65(GTX120030) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](#).