

KIAA0323 antibody [C2C3], C-term

Cat. No. GTX120182

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 75 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.59 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human KIAA0323. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx120182 IHC-P Image

KIAA0323 antibody [C2C3], C-term detects KHNYN protein at cytosol and nucleus on U251 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded U251 xenograft.

KIAA0323 antibody [C2C3], C-term (GTx120182) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




GTx120182 WB Image

Sample (30 ug of whole cell lysate)

A: HCT116


7.5% SDS PAGE

GTx120182 diluted at 1:5000



GTx120182 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using KIAA0323(GTx120182) antibody at 1:500 dilution.



GTx120182 WB Image

Sample (50 ug of whole cell lysate)

A: Mouse brain

7.5% SDS PAGE

GTx120182 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).