

NRG2 antibody

Cat. No. GTX12019

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5µg/ml
Not tested in other applications.	

Calculated MW 91.68 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1% Na ₂ HPO ₄ , 0.45% NaCl, 2.5% BSA
Preservative	0.025% Thimerosal, 0.025% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	500 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human NRG2(289-306aa NGRKNSRLQFNKVKVEDA), identical to the related rat sequence, and different from the related mouse sequence by one amino acid.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX12019 WB Image

WB analysis of various samples using GTX12019 NRG2 antibody.

Lane 1 : rat liver tissue lysate at 50ug

Lane 2 : COLO320 whole cell lysate at 40ug

Lane 3 : SMMC whole cell lysate at 40ug

Lane 4 : SW620 whole cell lysate at 40ug

Dilution : 0.5 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).