

WDYHV1 antibody

Cat. No. GTX120298

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.58 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human WDYHV1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX120298 IHC-P Image

WDYHV1 antibody detects WDYHV1 protein at nucleus in human breast cancer by immunohistochemical analysis.

Sample: Paraffin-embedded human breast cancer.

WDYHV1 antibody (GTX120298) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX120298 ICC/IF Image

WDYHV1 antibody detects WDYHV1 protein at cytoplasm and nucleus by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: WDYHV1 protein stained by WDYHV1 antibody (GTX120298) diluted at 1:200.

Blue: Hoechst 33342 staining.


GTX120298 IHC-P Image

WDYHV1 antibody detects WDYHV1 protein at cytosol and nucleus on mouse lung by immunohistochemical analysis.

Sample: Paraffin-embedded mouse lung.

WDYHV1 antibody (GTX120298) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](#).