

Nop25 antibody

Cat. No. GTX120303

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Nop25. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx120303 IHC-P Image

Nop25 antibody detects Nop25 protein at nucleous on mouse intestine by immunohistochemical analysis.

Sample: Paraffin-embedded mouse intestine.

Nop25 antibody (GTx120303) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




GTx120303 WB Image

Sample (30 ug of whole cell lysate)

A: Raji

12% SDS PAGE

GTx120303 diluted at 1:1000




GTx120303 IHC-P Image

Nop25 antibody detects NOL12 protein at nucleous on H1299 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded H1299 xenograft.

Nop25 antibody (GTx120303) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx120303 WB Image

Sample (50 ug of whole cell lysate)

A: mouse brain

12% SDS PAGE

GTx120303 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).