

LEF1 antibody [1C3.1D10]

Cat. No. GTX12033

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB
Reactivity	Human

Package

50 µg

Applications

Application Note

WB: Use at a concentration of 1 - 5 µg/ml. Predicted MW 55-57 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	44 kDa. (Note)
Product Note	Recognizes all human LEF1 isoforms using Western blotting. The antibody does not cross-react with TCF1 or TCF4.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide: TGIPHPAIVT PQVKQEHPT DSDLMHVKPQ HEQRKEQEPK RP, corresponding to amino acids 256-297 of Human LEF1.
Purification	Protein G purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).