

## RAB11FIP4 antibody

Cat. No. GTX120604

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | IHC-P      |
| Reactivity   | Mouse      |

Package  
100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.


## Properties

|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | PBS, 20% Glycerol                                                                                                                                                                                                                                                                                                                           |
| Preservative  | 0.01% Thimerosal                                                                                                                                                                                                                                                                                                                            |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                    |
| Immunogen     | Recombinant protein encompassing a sequence within the center region of human RAB11FIP4. The exact sequence is proprietary.                                                                                                                                                                                                                 |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX120604 IHC-P Image

RAB11FIP4 antibody detects RAB11FIP4 protein on mouse testis by immunohistochemical analysis.

Sample: Paraffin-embedded mouse testis.

RAB11FIP4 antibody (GTX120604) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).