

## Desmoglein 2 antibody [6D8]

## Cat. No. GTX12078

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	WB, ICC/IF, IP
<b>Reactivity</b>	Human

Package  
50 µg

## Applications

## Application Note

IF: Use at a dilution of 1/10. IP: Use at an assay dependent dilution. WB: Use at a dilution of 1/10. Predicted molecular weight: 165 kDa. Optimal dilutions/concentrations should be determined by the end user.

**Calculated MW** 122 kDa. ([Note](#))

**Product Note** GTX12078 reacts with Desmoglein 2. Does not cross-react with Desmoglein 1 or Desmoglein 3.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.1% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Unfortunately, this information is considered to be commercially sensitive
<b>Purification</b>	Protein G purified
<b>Conjugation</b>	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).