

DYNC1LI1 antibody

Cat. No. GTX120806

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:10000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human DYNC1LI1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx120806 WB Image


DYNC1LI1 antibody detects DYNC1LI1 protein by Western blot analysis.

A. 30 µg 293T whole cell lysate/extract

B. 30 µg HepG2 whole cell lysate/extract

7.5 % SDS-PAGE

DYNC1LI1 antibody (GTx120806) dilution: 1:1000



GTx120806 WB Image

DYNC1LI1 antibody detects DYNC1LI1 protein by western blot analysis.

A. 50 µg mouse brain lysate/extract

7.5 % SDS-PAGE

DYNC1LI1 antibody (GTx120806) dilution: 1:5000


GTx120806 ICC/IF Image

DYNC1LI1 antibody detects DYNC1LI1 protein at cytoplasm by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde/PBS for 15 min.

Green: DYNC1LI1 protein stained by DYNC1LI1 antibody (GTx120806) diluted at 1:500.

Blue: Hoechst 33342 staining.



For full product information, images and publications, please visit our [website](https://www.genetex.com).