

# RAB11-FIP5 antibody [N2C1], Internal

**Cat. No. GTX120985**

|                     |              |
|---------------------|--------------|
| <b>Host</b>         | Rabbit       |
| <b>Clonality</b>    | Polyclonal   |
| <b>Isotype</b>      | IgG          |
| <b>Applications</b> | WB, IHC-P    |
| <b>Reactivity</b>   | Human, Mouse |

**Package**  
100 µl, 25 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:3000         |
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.

**Calculated MW** 70 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | 0.1M Tris, 0.1M Glycine, 20% Glycerol                                                                                                                                                                                      |
| <b>Preservative</b>  | 0.01% Thimerosal                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.75 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                |
| <b>Immunogen</b>     | Recombinant protein encompassing a sequence within the center region of human RAB11-FIP5. The exact sequence is proprietary.                                                                                               |
| <b>Purification</b>  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


## DATA IMAGES



### GTx120985 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human hepatoma, using RAB11-FIP5(GTx120985) antibody at 1:500 dilution.


Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



### GTx120985 IHC-P Image

Immunohistochemical analysis of paraffin-embedded C2C12 xenograft, using RAB11-FIP5(GTx120985) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



### GTx120985 WB Image

Sample (30 µg of whole cell lysate)

A: 293T

B: A431

7.5% SDS PAGE

GTx120985 diluted at 1:1000

The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).