

TLR2 antibody

Cat. No. GTX12118

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 3 µg/ml
IHC-P	5 µg/ml

Not tested in other applications.

Calculated MW 90 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	This antibody was developed against a mixture of synthetic peptides containing amino acids 180-196, 353-370, and 473-489 of human TLR2 (NP_003255).
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX12118 WB Image

WB analysis of overexpressing TLR2 cell lysate using GTX12118 TLR2 antibody.



GTX12118 IHC-P Image

IHC-P analysis of mouse pancreas tissue using GTX12118 TLR2 antibody.

TOP : isotype control

Bottom : primary antibody

Dilution : 5 µg/ml



For full product information, images and publications, please visit our [website](#).