

PDIA6 antibody [N1N3]

Cat. No. GTX121275

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human PDIA6. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx121275 IHC-P Image

PDIA6 antibody [N1N3] detects PDIA6 protein at cytoplasm in mouse liver by immunohistochemical analysis.

Sample: Paraffin-embedded mouse liver.

PDIA6 antibody [N1N3] (GTx121275) diluted at 1:500.


Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx121275 IP Image

Immunoprecipitation of PDIA6 protein from 293T whole cell extracts using 5 µg of PDIA6 antibody [N1N3] (GTx121275).

Western blot analysis was performed using PDIA6 antibody [N1N3] (GTx121275).

EasyBlot anti-Rabbit IgG (GTx221666-01) was used as a secondary reagent.



GTx121275 IHC-P Image

PDIA6 antibody [N1N3] detects PDIA6 protein at cytoplasm in rat liver by immunohistochemical analysis.

Sample: Paraffin-embedded rat liver.

PDIA6 antibody [N1N3] (GTx121275) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx121275 WB Image

Sample (30 ug of whole cell lysate)

A: Raji

10% SDS PAGE

GTx121275 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).