

CD168 / RHAMM antibody [N1N2], N-term

Cat. No. GTX121502

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP, IHC
Reactivity	Human, Mouse, Rat

References (7)

★★★★☆ Review (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent
IHC-P	1:100-1:1000
IP	1:100-1:500
IHC	Assay dependent

Not tested in other applications.

Calculated MW 84 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human CD168 / RHAMM. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx121502 WB Image

Sample (20 ug of whole cell lysate)

A: human ESC

7.5% SDS PAGE

GTx121502 diluted at 1:1000



GTx121502 IHC-P Image


RHAMM antibody [N1N2], N-term detects RHAMM protein at cytosol on human breast carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human breast carcinoma.

RHAMM antibody [N1N2], N-term (GTx121502) dilution: 1:250.

Scale bar = 100 μ m.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx121502 IP Image

Immunoprecipitation of RHAMM protein from 293T whole cell extracts using 5 μ g of RHAMM antibody [N1N2], N-term (GTx121502).

Western blot analysis was performed using RHAMM antibody [N1N2], N-term (GTx121502).

EasyBlot anti-Rabbit IgG (GTx221666-01) was used as a secondary reagent.



For full product information, images and publications, please visit our [website](https://www.genetex.com).