

FSTL3 antibody

Cat. No. GTX12198

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5µg/ml
Not tested in other applications.	

Calculated MW 27.66 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1% Na ₂ HPO ₄ , 0.45% NaCl, 2.5% BSA
Preservative	0.025% Thimerosal, 0.025% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	500 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human FSTL3(145-159aa ECELRAARCRGHPDL), different from the related rat and mouse sequences by one amino acid.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX12198 WB Image

WB analysis of various samples using GTX12198 FSTL3 antibody.

Lane 1 : human U2OS whole cell lysates

Lane 2 : human HeLa whole cell lysates

Lane 3 : human A549 whole cell lysates


Lane 4 : human PC-3 whole cell lysates

Lane 5 : human HepG2 whole cell lysates

Lane 6 : human K562 whole cell lysates

Dilution : 0.5 µg/mL

Loading : 50µg of sample under reducing conditions


GTX12198 WB Image

WB analysis of various samples using GTX12198 FSTL3 antibody.

Lane 1 : rat liver tissue lysates

Lane 2 : rat testis tissue lysates

Lane 3 : rat NRK whole cell lysates

Lane 4 : mouse liver tissue lysates

Lane 5 : mouse testis tissue lysates

Lane 6 : mouse lung tissue lysates

Dilution : 0.5 µg/mL

Loading : 50µg of sample under reducing conditions



For full product information, images and publications, please visit our [website](https://www.genetex.com).