

HIBADH antibody

Cat. No. GTX122994

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 35 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.01 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human HIBADH. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX122994 WB Image

Sample (30 ug of whole cell lysate)

A: A549

10% SDS PAGE


GTX122994 diluted at 1:10000



GTX122994 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human breast cancer, using HIBADH(GTX122994) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX122994 WB Image


HIBADH antibody detects HIBADH protein by Western blot analysis.

A. 30 ug Neuro2A whole cell lysate/extract

B. 30 ug C8D30 whole cell lysate/extract

10 % SDS-PAGE

HIBADH antibody (GTX122994) dilution: 1:5000



GTX122994 WB Image

HIBADH antibody detects HIBADH protein by Western blot analysis.

A. 30 µg PC-12 whole cell lysate/extract

10 % SDS-PAGE

HIBADH antibody (GTX122994) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).