

CERCAM antibody

Cat. No. GTX123293

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 68 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human CERCAM. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX123293 WB Image

CERCAM antibody detects CERCAM protein by Western blot analysis.


A. 30 µg GL261 whole cell lysate/extract

B. 30 µg C8D30 whole cell lysate/extract

C. 30 µg NIH-3T3 whole cell lysate/extract


7.5 % SDS-PAGE

CERCAM antibody (GTX123293) dilution: 1:500


GTX123293 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human hepatoma, using CERCAM(GTX123293) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX123293 WB Image

CERCAM antibody detects CERCAM protein by Western blot analysis.

A. 50 µg mouse brain lysate/extract

7.5 % SDS-PAGE

CERCAM antibody (GTX123293) dilution: 1:500


GTX123293 IHC-P Image

CERCAM antibody detects CERCAM protein at cytoplasm in rat liver by immunohistochemical analysis.

Sample: Paraffin-embedded rat liver.

CERCAM antibody (GTX123293) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).