

14-3-3 beta antibody [4E1]

Cat. No. GTX12341

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IP
Reactivity	Mouse, Rat

Package
50 µg

Applications

Application Note

IP: Use at a concentration of 5 - 10 µg/ml. WB: Use at a concentration of 1 - 5 µg/ml. Detects a band of approximately 29 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% Sucrose
Preservative	1mg/ml Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein (Rat).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX12341 WB Image

Western blot analysis of mouse brain tissue extract, probed with GTX12341

Lane 1: Mouse IgG

Lane 2: Anti-14-3-3β



For full product information, images and publications, please visit our [website](#).