

HSP25 antibody

Cat. No. GTX12351

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Mouse, Rat, Hamster, Guinea pig, Pig

References (2)

Package

50 µl

Applications

Application Note

For IP: Use at an assay dependent dilution. For WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant mouse Hsp25
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx12351 WB Image

Western Blot analysis of HSP25, pAb.

Lane 1: MW Marker

Lane 2: PC-12 Cell Lysate (Heat Shock)

Lane 3: 3T3 Cell Lysate (Heat Shock)

Lane 4: HSP25 (mouse), (recombinant) (GTx12351)



For full product information, images and publications, please visit our [website](https://www.genetex.com).