

alpha Tubulin antibody

Cat. No. GTX124303

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Wm, IHC
Reactivity	Zebrafish

References (4) Package 100 μΙ, 25 μΙ

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
IHC-Wm	1:100-1:500
IHC	Assay dependent

Not tested in other applications.

Properties	
Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.24 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen used to generate this antibody corresponds to zebrafish alpha Tubulin
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 02 Page 1 of 2




DATA IMAGES



GTX124303 IHC-Wm Image

alpha Tubulin antibody detects alpha Tubulin protein on zebrafish by whole mount immunohistochemical analysis.

Sample: 2 days-post-fertilization zebrafish embryo. alpha Tubulin antibody (GTX124303) dilution: 1:200.



GTX124303 IHC-P Image


Immunohistochemical analysis of paraffin-embedded zebrafish mature sperms, using alpha Tubulin (GTX124303) antibody at 1:300 dilution.



GTX124303 IHC-Wm Image

alpha Tubulin antibody detects alpha Tubulin protein on zebrafish by whole mount immunohistochemical analysis.

Sample: 2 days-post-fertilization zebrafish embryo. alpha Tubulin antibody (GTX124303) dilution: 1:200.



GTX124303 WB Image

Sample (10 ug of whole cell lysate) A: Adult zebrafish 10% SDS PAGE GTX124303 diluted at 1:1000



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 02 Page 2 of 2