

FGF1 antibody (Biotin)

Cat. No. GTX12474

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

Package

25 µg

Applications

Application Note

ELISA: Use at a concentration of 0.15 - 0.3 µg/ml. Sensitivity = 0.2 ng/wellWB: Use at a concentration of 0.1 - 0.2 µg/ml. Predicted molecular weight: 18.6 kDa. Sensitivity = 1.5-3.0 ng/lane Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 17 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein (Human)
Purification	Immunogen affinity purified Antiserum was affinity purified using immobilized recombinant Human FGF acidic.
Conjugation	Biotin
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).