

MCHR2 antibody

Cat. No. GTX12674

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	10 - 20 µg/ml
IHC-P	21 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 19 amino acid peptide from C-terminus of human MCHR2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX12674 IHC-P Image


IHC-P analysis of human brain, neurons tissue using GTX12674 MCHR2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX12674 IHC-P Image

IHC-P analysis of brain, amygdala, neurons and cell processes tissue using GTX12674 MCHR2 antibody.


GTX12674 IHC-P Image

IHC-P analysis of brain, caudate, neurons and glia tissue using GTX12674 MCHR2 antibody.


GTX12674 ICC/IF Image

ICC/IF analysis of HEK293 transfected with MCHR2 using GTX12674 MCHR2 antibody.



For full product information, images and publications, please visit our [website](#).