

MYT1 antibody

Cat. No. GTX128188

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 122 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.96 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human MYT1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx128188 WB Image


Myt1 antibody detects Myt1 protein by western blot analysis.

A. 30 μ g IMR32 whole cell extract

B. 30 μ g IMR32 nuclear extract

5 % SDS-PAGE

Myt1 antibody (GTx128188) dilution: 1:500



GTx128188 WB Image

Myt1 antibody detects Myt1 protein by western blot analysis.

A. 30 μ g SK-N-SH nuclear extract

5 % SDS-PAGE


Myt1 antibody (GTx128188) dilution: 1:500


GTx128188 IHC-P Image

Myt1 antibody detects Myt1 protein at nucleus on mouse skin by immunohistochemical analysis.

Sample: Paraffin-embedded mouse middle brain.

Myt1 antibody (GTx128188) dilution: 1:500.


GTx128188 WB Image

Myt1 antibody detects MYT1 protein by Western blot analysis.

A. 30 μ g U87-MG whole cell lysate/extract

5 % SDS-PAGE

Myt1 antibody (GTx128188) dilution: 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).