

USP3 antibody

Cat. No. GTX128238

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

[References \(6 \)](#)

[Package](#)

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:20000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 59 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human USP3. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX128238 WB Image

USP3 antibody detects USP3 protein by western blot analysis.


A. 30 μ g 293T whole cell lysate/extract

B. 30 μ g A431 whole cell lysate/extract

C. 30 μ g HeLa whole cell lysate/extract

7.5 % SDS-PAGE

USP3 antibody (GTX128238) dilution: 1:1000



GTX128238 ICC/IF Image

USP3 antibody detects USP3 protein at nucleus by immunofluorescent analysis.

Sample: K562 cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: USP3 protein stained by USP3 antibody (GTX128238) diluted at 1:500.

Blue: Hoechst 33342 staining.


GTX128238 WB Image


USP3 antibody detects USP3 protein by western blot analysis.

A. 30 μ g 293T whole cell lysate/extract

B. 30 μ g whole cell lysate/extract of human USP3-transfected 293T cells

7.5 % SDS-PAGE

USP3 antibody (GTX128238) dilution: 1:20000


GTX128238 WB Image

The data was published in the journal Oncotarget in 2017. [PMID: 28881649](https://pubmed.ncbi.nlm.nih.gov/28881649/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).