

MMP24 antibody

Cat. No. GTX128246

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000

Not tested in other applications.

Calculated MW 73 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the Extracellular domain of human MMP24. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx128246 WB Image


MMP24 antibody detects MMP24 protein by Western blot analysis.

A. 30 µg 293T whole cell lysate/extract

B. 30 µg whole cell lysate/extract of human MMP24-transfected 293T cells

7.5 % SDS-PAGE

MMP24 antibody (GTx128246) dilution: 1:5000




GTx128246 WB Image

MMP24 antibody detects MMP24 protein by Western blot analysis.

A. 50 µg mouse brain lysate/extract

7.5 % SDS-PAGE

MMP24 antibody (GTx128246) dilution: 1:1000



GTx128246 WB Image

MMP24 antibody detects MMP24 protein by western blot analysis.

A. 30 µg SK-N-SH whole cell lysate/extract

B. 30 µg IMR32 whole cell lysate/extract

7.5 % SDS-PAGE

MMP24 antibody (GTx128246) dilution: 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).