

## Malondialdehyde antibody

Cat. No. GTX12835

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA, Immunoblotting, IHC
Reactivity	Species independent

References ( 1 )

Package

50 µl

## PRODUCT

Summary	Compounds	Affinity compared to Malondialdehyde-BSA
	Malondialdehyde-BSA	1
	4-hydroxynonenal-BSA	> 100,000 Times
	Acrolein-BSA 1/> 100,000	> 100,000 Times
	Formaldehyde-BSA	> 100,000 Times
	Azelaic acid-BSA	> 100,000 Times

## Applications

## Application Note

Immunohistochemistry, Immunoblot. Optimal dilutions should be determined by each laboratory for each application.

## Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Purification	Antiserum previously preabsorbed on protein carriers, and purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2025 / 05 / 10 Page 1 of 2

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).