

Cullin 4b antibody

Cat. No. GTX129298

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 104 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.75 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Cullin 4b. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx129298 IHC-P Image

Cullin 4b antibody detects Cullin 4b protein at nucleus on human esophagus carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human esophagus carcinoma.

Cullin 4b antibody (GTx129298) dilution: 1:1000.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx129298 WB Image

Cullin 4b antibody detects Cullin 4b protein by western blot analysis.

A. 30 µg Neuro2A whole cell lysate/extract

B. 30 µg C8D30 whole cell lysate/extract


C. 30 µg NIH-3T3 whole cell lysate/extract

D. 30 µg Raw 264.7 whole cell lysate/extract

E. 30 µg C2C12 whole cell lysate/extract

5 % SDS-PAGE

Cullin 4b antibody (GTx129298) dilution: 1:1000




GTx129298 WB Image

Cullin 4b antibody detects Cullin 4b protein by western blot analysis.

A. 30 µg A375 whole cell lysate/extract

5 % SDS-PAGE

Cullin 4b antibody (GTx129298) dilution: 1:1000



GTx129298 WB Image

Cullin 4b antibody detects Cullin 4b protein by western blot analysis.

A. 30 µg PC-12 whole cell lysate/extract

B. 30 µg Rat2 whole cell lysate/extract

5 % SDS-PAGE

Cullin 4b antibody (GTx129298) dilution: 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).