

TSLP antibody

Cat. No. GTX129591

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
ELISA	1:-1:10000

Not tested in other applications.

Calculated MW 18 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human TSLP. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx129591 ELISA Image

An ELISA plate is coated with 50 μ L of TSLP recombinant protein at concentration of 2 μ g/mL. The coated protein is detected with anti-TSLP antibody (GTx129591) at concentration ranged from 31.25 to 500 ng/mL.



GTx129591 WB Image

TSLP antibody detects TSLP protein by western blot analysis. Human tissue extracts (30 μ g) was separated by 15 % SDS-PAGE, and the membrane was blotted with TSLP antibody (GTx129591) diluted by 1:500.



GTx129591 IHC-P Image

TSLP antibody detects TSLP protein at cytoplasm on human colon carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human colon carcinoma.

TSLP antibody (GTx129591) diluted at 1:100.



For full product information, images and publications, please visit our [website](https://www.genetex.com).