

## MED13 antibody

Cat. No. GTX129674

|              |              |
|--------------|--------------|
| Host         | Rabbit       |
| Clonality    | Polyclonal   |
| Isotype      | IgG          |
| Applications | WB           |
| Reactivity   | Human, Mouse |

References ( 4 )

Package

100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution |
|-----------------------------------|----------------------|
| WB                                | 1:500-1:3000         |
| Not tested in other applications. |                      |

**Calculated MW** 239 kDa. ( [Note](#) )**Product Note** Knockdown/Knockout validation was supported by references (PMID:36220895)


## Properties

|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | PBS, 20% Glycerol                                                                                                                                                                                                                                                                                                                           |
| Preservative  | 0.025% ProClin 300                                                                                                                                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                    |
| Immunogen     | Recombinant protein encompassing a sequence within the center region of human MED13. The exact sequence is proprietary.                                                                                                                                                                                                                     |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




**GTX129674 WB Image**

Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with MED13 antibody (GTX129674) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



**GTX129674 WB Image**

The data was published in the 2018 in Elife. [PMID: 29809150](https://doi.org/10.7554/eLife.29809)



**GTX129674 WB Image**

The data was published in the 2018 in Elife. [PMID: 29809150](https://doi.org/10.7554/eLife.29809)



For full product information, images and publications, please visit our [website](https://www.genetex.com).