

Rho GDI alpha antibody

Cat. No. GTX13062

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse, Bovine

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-1.0 µg/ml
IP	4 µg

Not tested in other applications.

Calculated MW	23 kDa. (Note)
Product Note	The antibody recognizes Rho-GDI and cross-reacts with D4-GDI.

Properties

Form	Liquid
Buffer	100mM Tris-Glycine, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	~1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GST fusion protein corresponding to full-length human Rho-GDI.
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).