

BTF antibody

Cat. No. GTX130772

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 106 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.08 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human BTF . The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx130772 WB Image

K562 whole cell and nuclear extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with BTF antibody (GTx130772) diluted at 1:1000.



GTx130772 ICC/IF Image


BTF antibody detects BTF protein at nucleus by immunofluorescent analysis.
Sample: MCF7 cells were fixed in 4% paraformaldehyde at RT for 15 min.
Green: BTF protein stained by BTF antibody (GTx130772) diluted at 1:500.
Blue: Hoechst 33342 staining.
Scale bar = 10 µm.



GTx130772 IHC-P Image

BTF antibody detects BTF protein at nucleus in rat brain by immunohistochemical analysis.
Sample: Paraffin-embedded rat brain.
BTF antibody (GTx130772) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx130772 IHC-P Image

BTF antibody detects BTF protein at nucleus in mouse cervix by immunohistochemical analysis.
Sample: Paraffin-embedded mouse cervix.
BTF antibody (GTx130772) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).