

# PTPRA antibody

## Cat. No. GTX13095

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rat, Zebrafish, Rabbit, Bovine, Dog, Hamster, Chicken, Pig, Monkey, Bat, Horse, Xenopus

Package 25 μg

## Applications

#### **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	4 μg/ml

Not tested in other applications.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 18 amino acid peptide from cytoplasmic domain of human PTPRA.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 2



### DATA IMAGES



#### GTX13095 IHC-P Image

IHC-P analysis of human breast, epithelium tissue using GTX13095 PTPRA antibody.

Antigen retrieval: Heat-induced antigen retrieval

Dilution: 4 μg/ml



#### GTX13095 IHC-P Image

 $IHC-P\ analysis\ of\ human\ skeletal\ muscle\ tissue\ using\ GTX13095\ PTPRA\ antibody.$ 

Antigen retrieval: Heat-induced antigen retrieval

Dilution :  $4 \, \mu g/ml$ 

For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 2 of 2