

proBNP antibody [5B6]

Cat. No. GTX13111

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package

200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	15 kDa. (Note)
Product Note	NT-proBNP and proBNP

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic N-terminal pro brain natriuretic peptide (NT-proBNP), corresponding to amino acid residues 1-12. The peptides used for immunization were conjugated with carrier protein.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX13111 WB Image

Detection of human recombinant NT-proBNP expressed in *E. coli* in Western blotting by different monoclonal antibodies after Tricine-SDS gel electrophoresis. Antigen: recombinant NT-proBNP expressed in *E. coli*, 2.5 ug; Lane 1: proBNP antibody [5B6] (GTX13111); Lane 2: proBNP antibody [15F11] (GTX13115); Lane 3: proBNP antibody [11D1] (GTX13059); Lane 4: proBNP antibody [24E11] (GTX13123)



GTX13111 WB Image

Detection of human recombinant proBNP expressed mammalian cells in Western blotting by different monoclonal antibodies after Tricine-SDS gel electrophoresis. Antigen: recombinant proBNP expressed in mammalian cells, 0.15 ug; Lane 1: proBNP antibody [5B6] (GTX13111); Lane 2: proBNP antibody [15F11] (GTX13115); Lane 3: proBNP antibody [3G12] (GTX14712); Lane 4: proBNP antibody [11D1] (GTX13059); Lane 5: proBNP antibody [24E11] (GTX13123); Both proBNP antibody [11D1] (GTX13059) and proBNP antibody [24E11] (GTX13123) do not recognize glycosylated proBNP expressed in mammalian cells.



For full product information, images and publications, please visit our [website](https://www.genetex.com).