

Adenosine A2b Receptor antibody

Cat. No. GTX132217

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.53 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human Adenosine A2b Receptor. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX132217 IHC-P Image

Adenosine A2b Receptor antibody detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody (GTX132217) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX132217 IHC-P Image

Adenosine A2b Receptor antibody detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse cerebellum.

Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody (GTX132217) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX132217 IHC-P Image

Adenosine A2b Receptor antibody detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat cerebellum.

Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody (GTX132217) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX132217 IHC-P Image

Adenosine A2b Receptor antibody detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat brain.

Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody (GTX132217) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).