

Bassoon antibody [SAP7F407]

Cat. No. GTX13249

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-Fr, IP, IHC
Reactivity	Human, Mouse, Rat

References (9)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	Assay dependent
IHC-Fr	Assay dependent
IP	2.5µg
IHC	1:400

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol, 0.1mM PMSF
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant rat bassoon fragment (residues 738-1035) expressed as a fusion protein in E.coli.
Purification	Protein G purified This antibody has been affinity purified.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx13249 IHC Image

Immunohistochemistry analysis of bassoon localization in the rat brain cerebellar molecular layer using Bassoon, mAb (SAP7F407)



GTx13249 WB Image

Western blot analysis of rat brain tissue lysate using Bassoon antibody [SAP7F407].



For full product information, images and publications, please visit our [website](#).