

Valproic Acid antibody

Cat. No. GTX13257

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Species independent

Package

50 µl

PRODUCT

Summary

Compounds	Affinity compared to Valproic acid-BSA
Valproic acid-BSA	1
Methylpentanoic acid-BSA	33 Times
Caproic acid-BSA	>10,000 Times
Valeric acid-BSA	>10,000 Times
Butyric acid-BSA	>10,000 Times
Free Valproic acid	>10,000 Times

Applications

Application Note

Optimal dilutions should be determined by each laboratory for each application

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	No preservatives



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic Valproic acid conjugated to Protein Carriers (PC).
Purification	previously preabsorbed on protein carriers, and purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).