

Morganella morganii antibody

Cat. No. GTX13265

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Morganella morganii

Package

50 µl

PRODUCT

Summary

Compounds	Affinity compared to Morganella morganii
Morganella morganii	1
Hafnia alvei	50 Times
Pseudomonas aeruginosa	50 Times
Pseudomonas putida	50 Times
Proteus mirabilis	50 Times
Citrobacter koseri	50 Times
Klebsiella pneumoniae	50 Times

Applications

Application Note

Optimal dilutions should be determined by each laboratory for each application

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



For full product information, images and publications, please visit our [website](#).

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Morganella morganii total proteins
Purification	Ammonium Sulphate purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).