

## GPRC5A antibody

Cat. No. GTX13274

|              |               |
|--------------|---------------|
| Host         | Rabbit        |
| Clonality    | Polyclonal    |
| Isotype      | IgG           |
| Applications | ICC/IF, IHC-P |
| Reactivity   | Human         |

## Package

25 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ICC/IF             | Assay dependent      |
| IHC-P              | 8 - 38 µg/ml         |

Not tested in other applications.

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.1% Sodium azide                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Synthetic 16 amino acid peptide from N-terminal extracellular domain of human GPRC5A / RAI3.                                                                                                                               |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX13274 ICC/IF Image

ICC/IF analysis of HEK293 transfected with GPRC5A using GTX13274 GPRC5A antibody.



### GTX13274 IHC-P Image

IHC-P analysis of human breast tissue using GTX13274 GPRC5A antibody.

Antigen retrieval : Heat-induced antigen retrieval



For full product information, images and publications, please visit our [website](#).