

BAI2 antibody

Cat. No. GTX13306

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	IHC-P
Reactivity	Human, Mouse, Rat, Dog, Hamster, Pig, Monkey, Horse

Package

25 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	2 - 3 µg/ml
Not tested in other applications.	

Calculated MW 173 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 16 amino acid peptide from N-terminal extracellular domain of human BAI2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX13306 IHC-P Image


IHC-P analysis of human ovary, carcinoma tissue using GTX13306 BAI2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13306 IHC-P Image


IHC-P analysis of human breast, carcinoma tissue using GTX13306 BAI2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13306 IHC-P Image

IHC-P analysis of human lung, small cell carcinoma tissue using GTX13306 BAI2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13306 IHC-P Image

IHC-P analysis of human brain, neurons and glia tissue using GTX13306 BAI2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13306 IHC-P Image

IHC-P analysis of brain, glioblastoma tissue using GTX13306 BAI2 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).