

Leukotriene B4 Receptor 2 antibody

Cat. No. GTX13348

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr, FCM
Reactivity	Human, Monkey

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	26 µg/ml
IHC-Fr	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 16 amino acid peptide from 3rd cytoplasmic domain of human LTB4R2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX13348 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX13348 Leukotriene B4 Receptor 2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13348 IHC-P Image

IHC-P analysis of human lymph node, hodgkins lymphoma tissue using GTX13348 Leukotriene B4 Receptor 2 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX13348 IHC-P Image

IHC-P analysis of human colon, carcinoma and leukocytes tissue using GTX13348 Leukotriene B4 Receptor 2 antibody.
Antigen retrieval : Heat-induced antigen retrieval



For full product information, images and publications, please visit our [website](https://www.genetex.com).