

Sarcosine antibody

Cat. No. GTX13479

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Species independent

Package
50 µl

PRODUCT

Summary

Compounds	Affinity compared to Sarcosine-BSA
Sarcosine-BSA	1
Glycine-BSA	>50,000 Times

Applications

Application Note

Optimal dilutions should be determined by each laboratory for each application.

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic Sarcosine conjugated to protein carriers
Purification	previously preabsorbed on protein carriers, and ammonium sulphate purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).