

MATH1 (188-200) peptide (BSA)

Cat. No. GTX13519

Species Human

Package
500 µg

Properties

Form Liquid

Buffer PBS, 40% Glycerol

Preservative No preservatives

Storage Upon receipt - Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Synthetic peptide: SFNNDKCLKSKYET, corresponding to amino acids 188-200 of Human MATH1.

Purification Purified

Conjugation BSA

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).